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- Gallo RC, Salahuddin SZ, Popovic M, et al. Frequent detection and isolation of cytopathic retroviruses (HTLV III) from patients with AIDS and at risk for AIDS. Science 1984; 224:500-3.
- Curran JW, Morgan WM, Hardy AM, et al. The epidemiology of AIDS: current status and future prospects. Science 1985; 222:352-7.
- Clavel F, Guetard D, Brun-Vezinet F, et al. Isolation of a new human retrovirus from West African patients with AIDS. Science 1986; 233:343-6.
- Rosenberg ZF and Fauci AS. The immunopathogenesis of HIV infection. Frank J. Dixon, Editor. Advances in Immunology, 1989; 47:377-431.
- Schochetman G, and George JR., Editors. AIDS Testing Methodology and Management Issues. Springer-Verlag, New York City, N.Y., 1992.
- Farr J, Perry KR, Mortimer PP. Sensitive assays for viral antibodies in saliva: an alternative to tests on serum. Lancet 1987; 335:72-75.
- Holmstrom P, Svanen S, Laine P, Valle SL, Suni J. HIV antibodies in whole saliva detected by ELISA and western blot assays. Journal of Medical Virology 1990; 30(4):245-248.
- Holm-Hansen C, Constantine NT, Haukenes G. Detection of antibodies to HIV in homologous sets of plasma, urine and oral mucosal transudate samples using rapid assays in Tanzania. Clinical and Diagnostic Virology 1993; 1(3):207-214.
- Granade TC, Phillips SK, Parash B, Pau CP and George JR. Oral fluid as a specimen for detection and confirmation of antibodies to human immunodeficiency virus type 1. Clinical and Diagnostic Laboratory Immunology 1995; 2(4):395-399.
- Gottfried, TD, Mink RW and Phanuphak P. Calypte Aware HIV-1/2 OMT antibody test using oral fluid: special challenges of rapid HIV testing in the developing world. Expert Rev. Mol. Diagn. 2006; 6(2):139-144.
- CDC. Recommendations for prevention of HIV transmission in health-care settings. MMWR 1987; 36(SU02):3S-18S.
- CDC. Perspectives in disease prevention and health promotion update: universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus, and other bloodborne pathogens in health-care settings. MMWR 1988; 37(24):377-388.
- CDC. Guidelines for prevention of transmission of human immunodeficiency virus and hepatitis B virus to healthcare and public-safety workers. MMWR 1989; 38(S-6):3-37.
- Sehulster LM, Hollinger FB, Dressman GR, and Melnick JL. Immunological and biophysical alteration of hepatitis B virus antigens by sodium hypochlorite disinfection. Appl Env Microbiol 1981; 42: 762-7.

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